

# Validation of Mother Stock 1 Panel Duplex-TaqMan qPCR Assay

## Background



Starting with clean, tested stock is a critical step in growing a successful crop. Among cannabis pathogens, hop latent viroid and Fusarium are by far the most economically damaging infections. Beginning a run with mother stock infected with either of these two pathogens can lead to catastrophic failure of a harvest and potentially facility-wide spread. This assay detects both of these pathogens simultaneously to allow economical screening of the most dangerous cannabis pathogens in a single sample.

Hop Latent Viroid (HLVd) is a small, circular, infectious agent. As opposed to viruses and other pathogens, HLVd lacks an outside protein and is only composed of genetic material (RNA). Viroids generally spread via mechanical transmission on unsterilized cutting tools and equipment. HLVd has been detected in most geographical locations around the world and identified in cannabis plants throughout the United States. Common symptoms of HLVd include stunted growth, brittle stems, leaf malformation and reduced flower mass. However, plants may initially appear asymptomatic or with subtle symptoms making detection by eye difficult.

Fusarium is a root-infecting fungal pathogen. Fusarium infection can result in both root rot and dampingoff (i.e. stunted growth with dark green or yellow curled leaves and root lesions, especially at root tips). In an indoor growing environment, introduction of diseased plant material such as cuttings, stock plants or seeds can result in inoculation and eventual spread of the fungus. Infected growing medium that has not been adequately sterilized can also be a source of infection. Once established, Fusarium can spread via unfiltered, recirculated water, in air, on tools/equipment and potentially by workers

## General Information

TUMI Genomics' Mother Stock Screening diagnostic assay is for the detection of infection in cannabis and hemp plant tissue. The assay uses a duplex TaqMan qPCR to detect the presence of Fusarium DNA/RNA sequences and hop latent viroid RNA sequences. The primer/probe sets target the pathogen genomes as well as an endogenous plant transcript as a positive control. Because this assay detects a target that is in all cannabis tissues, it provides confirmation of sample integrity, successful nucleic acid extraction, and activity of reaction components in every sample.



## Test Validation

**Technical validation:** Technical validation consists of experiments performed by TUMI Genomics to determine the technical (in lab) limits and characteristics of the assay.

**Limit of Detection (LoD):** A limit of detection study determines the sensitivity of a test. The limit of detection is defined as the lowest concentration of pathogen where the test can still detect 95% of the true positive samples. A limit of detection can help you understand how well a test can find low-level/early infections.

The LoD (lowest concentration of pathogen detected by this assay) was determined by testing serial dilutions (500 copies – 415copies) of genomic material spiked into negative cannabis plant extract and tested with ten replicates per concentration.

The smallest amount of both *Fusarium* and hop latent viroid reliably detected by the assay was determined to be ~ 30 genomic copies per reaction.

The table below shows the results of tests to determine the limit of detection (i.e. the smallest amount of pathogen this assay can identify). The limit of detection is highlighted in green. The limit of detection was confirmed for each individual species of *Pythium* as well.

Sample Concentration	Fraction Positive	Mean CT Plant Control	Standard Deviation	Mean CT HLVd (Target #1)	Standard Deviations	Mean CT HLVd (Target #2)	Standard Deviations
500 copies/reaction	10/10	28.72	+/- 1.04	26.81	+/- 0.13	27.68	+/- 0.32
250 copies/reaction	10/10	29.07	+/- 1.32	27.08	+/- 0.11	28.13	+/- 0.21
125 copies/reaction	10/10	29.56	+/- 1.57	28.38	+/- 0.16	29.61	+/- 0.23
31 copies/reaction	10/10	30.24	+/- 0.14	30.44	+/- 0.87	31.05	+/- 0.15
15 copies/reaction	10/10	29.65	+/- 0.07	31.41	+/- 0.29	32.06	+/- 0.71
7.5 copies/reaction	10/10	29.42	+/- 0.16	32.01	+/- 0.122	33.03	+/- 0.45
3.8 copies/reaction	3/10	29.23	+/- 0.15	32.71	+/- 1.53	33.93	+/- 0.73
0 copies/reaction	0/30	28.76	+/- 1.16	NaN	NaN	NaN	NaN

**Cross Reactivity:** A cross reactivity analysis determines the specificity of a test. Specificity means whether the test detects only the target pathogen versus giving a non-specific signal or throwing a positive result due to the presence of an unrelated pathogen. A specific test should be 100% specific for the target pathogen (i.e. hop latent viroid) and less than 80% specific for potential contaminating pathogens (like other viruses).

Cross-reactivity was evaluated using in silico analysis of primer sequences compared to genome sequences from microorganisms that are commonly found in the roots, stems, and leaves of Cannabis, hops, and hemp plants. Genomic sequences for the microorganisms were acquired from the National Center for Biotechnology Information database (NCBI) and the alignments were performed with the Basic local alignment search tool (BLAST). For each organism, percent cross-reactivity was determined by dividing the number of nucleotide sequences that matched with the organism by the total number of nucleotides in the *Pythium* primer/probe set.

Cross-reactivity is defined as greater than 80% similarity between the primer/probes set and any sequence present in the targeted microorganism. As expected, primer probe sequences matched sequences from the *Pythium* genome 100%. No cross-reactivity above 80% was found with any other tested microorganism.

For details from this analysis please reference the following two documents: Hop latent viroid triplex TaqMan qRT-PCR validation [Fusarium duplex TaqMan qPCR validation](#)

**Inclusivity:** An inclusivity analysis determines how well a test can detect different known sequence variants of a given pathogen. An acceptable inclusivity analysis should show that an assay has been designed in a way that can reasonably detect all, or the vast majority, of known variants of the target pathogen.

This analysis indicated that the primers and probes used in this assay reliably detect 98% of known hop latent viroid variants (152/155) with little to no reduction in amplification efficiency and will identify 100% of all known hop latent viroid variants in a yes/no diagnostic assay.

To ensure inclusivity of the assay, all *Fusarium* primer and probe sequences were designed specifically against reference genomes identified from infections in cannabis (reviewed in Gwinn et. al, 2022). All accessions identified from pathogenic cannabis tissue can be amplified using the primer sequences.

\*Gwinn, Kimberly et. al, Diseases of Cannabis sativa Caused by Diverse Fusarium Species. *Frontiers in Agronomy*, Jan 18, 2022. <https://doi.org/10.3389/fagro.2021.796062>

## **TUMI Genomics Lead Scientists**

**Tassa Saldi, PhD:** Dr. Saldi received her undergraduate and graduate degrees in molecular biology from the University of Colorado in Boulder and completed her post-doctoral studies at the Health Sciences Center, University of Colorado, Denver. Her graduate work explored the molecular mechanism underpinning Amyotrophic Lateral Sclerosis (ALS) and the role of double-stranded RNA accumulation and heterochromatin in pathogenesis.

Continuing her work on structured RNA during post-doctoral work, Dr. Saldi investigated the role of genome-wide nascent RNA secondary structure in co-transcriptional splicing, A-to-I RNA editing and transcription termination. Her work was supported by fellowships from the American Cancer Society and the RNA Biosciences Intuitive (RBI). Following her postdoc, Dr. Saldi directed the COVID-19 surveillance lab at CU, Boulder where she supervised a team of 8 scientists and designed and validated multiple PCR assays to detect SARS-CoV-2 in human saliva. She is a lead scientist and CSO of TUMI Genomics.

Her publications can be found here: <https://pubmed.ncbi.nlm.nih.gov/term=Tassa+Saldi&sort=date>

**Alfonso Garrido-Lecca, PhD:** Dr. Garrido-Lecca received an undergraduate degree in biology with a minor in chemistry from Texas A&M University. He pursued his PhD at the University of Colorado, Boulder in molecular biology. His graduate work focused on using the unique genetic organization of *C. elegans* to understand how genes are expressed and RNA transcripts processed. His postdoctoral work focused on the regulation of microRNAs in leukemia and was supported by a fellowship from the Linda Crnic Institute for Down Syndrome and the National Institute of Health T32 training grant. Dr. GarridoLecca is a lead scientist at TUMI Genomics and head of Research and Development.

His publications can be found here:  
<https://pubmed.ncbi.nlm.nih.gov/term=alfonso+garridolecca&sort=date>